

# ACCESSION+ BY *DATARAY*

- Conveniently prints Send Patient collection labels in the patient draw room, reducing the risk of mixing up accession labels.
- Automatically generates Immediate Collect and Lab Collect labels to eliminate handwritten labels, and to avert the transport of unlabeled tubes.
- Easily prints Room Collect labels in the patient care rooms for non-blood items such as throat swabs, etc.
- Quickly initiates an accession by allowing you to simply scan the barcode the new collections have.

## IN ACTION

- 1 • A provider enters the collection order into EHR.
- 2 • The sample is collected and a DataRay Room Collect Label is attached to the vial.\* The collection is taken to the lab.
- 3 • The lab technician scans the barcode on the Room Collect label to access the patient's collection order in the LIS.
- 4 • Accession+ label is printed and attached to the collection to begin the accession process using the analyzer / biomedical equipment.

\* see DPOCC flyer



Accession+ Label Layout

**DEMO, DATARAY ADULT**  
999999 DOB: 01/01/2001  
CH 9999 9 04/03/19 10:50 LOC: ER ORD#: 12345 Female  
PROV: Dr. John Cleese  
MARBLE/SST CMP



- Compatible with virtually every EHR including VistA, RPMS, Epic, Meditech, and Cerner

